

3M™ Ultimate 2000 Veterinary X-ray Film

*The highest over/under exposure tolerance
in veterinary medicine ...*



Finally a Forgiving Film!

3M *Innovation*

3M™ Ultimate 2000 Veterinary X-ray Film

The latitude you want, the resolution you need


At 3M Animal Care, we understand the challenge of taking radiographs quickly, easily and effectively. Whether you're evaluating an injured dog or doing pre-purchase exams in the sale barn, you need to know that your first film is the only one you'll need. And, given the high costs of running your own veterinary business, you simply can't afford the time and the expense of multiple retakes.

3M™ **Ultimate 2000** Veterinary X-ray film is the one brand of film your practice can't afford to be without. Whether you're using your portable X-ray machine in the barn or in the clinic, 3M™ **Ultimate 2000** minimizes the effects of patient motion to deliver a perfect radiograph every time. This state-of-the-art, green-sensitive, high speed, medium contrast double emulsion film utilizes an exclusive dye underlayer to deliver maximum detail for a precise evaluation with minimal exposure. No matter how many patients you see in a day, you'll find that 3M™ **Ultimate 2000** is so precise, so consistent and so reliable, you'll only need to take one film for each view -- and only one brand of film for your veterinary practice.

Depend on 3M™ Animal Care to reduce your X-ray costs

- 3M Satisfaction Guarantee
- 3M Technique guides
- 3M Technical assistance
- 3M Technician toll-free telephone support
- Cassette acquisition programs

Double emulsion construction delivers ultimate contrast with minimum exposure



FEATURES	ADVANTAGES	BENEFITS
<ul style="list-style-type: none"> • Double emulsion construction 	<ul style="list-style-type: none"> • Decreased exposure time 	<ul style="list-style-type: none"> • Minimizes effects of patient motion, reduces unnecessary radiation exposure, extends X-ray tube life, can be used with portable machines
<ul style="list-style-type: none"> • Zero crossover technology 	<ul style="list-style-type: none"> • Prevents light crossover, improves detail, decreased 'noise' • Provides wide range of over/under exposure tolerance 	<ul style="list-style-type: none"> • Increased diagnostic information • Fewer retakes • Reduce cost • Reduce staff & patient exposure
<ul style="list-style-type: none"> • Wide latitude film 	<ul style="list-style-type: none"> • Provides balanced image visibility of soft tissue and trabecular bone 	<ul style="list-style-type: none"> • One film for all views
<ul style="list-style-type: none"> • Safe for all processing conditions 	<ul style="list-style-type: none"> • Convenient 	<ul style="list-style-type: none"> • No special equipment required



Animal Care Products
 3M Center, Building 270-2N-03
 St. Paul, MN 55144-1000
 1-800-848-0829

