

3M™ Ultra Detail Plus Veterinary X-ray Film

The resolution you want ...



The reliability you need!

3M *Innovation*

3M™ Ultra Detail Plus Veterinary X-ray Film


Technology to meet the needs of the Veterinary Professional

At 3M™ Animal Care, we understand the challenge of taking radiographs quickly, easily and effectively. Whether you're evaluating an injured dog or doing a pre-purchase exam on a horse, you need to know that you'll get it right the first time and every time; no repeats required. Given the high cost of running a veterinary business, you simply can't afford the time and expense of multiple retakes.


3M™ **Ultra Detail Plus** Veterinary X-ray Film is the film you have been looking for. 3M™ **Ultra Detail Plus** is a state-of-the-art, green-sensitive, high speed, high contrast, double emulsion film that utilizes an exclusive dye underlayer to deliver maximum detail at a minimal exposure. Whether you are using a portable X-ray machine in a barn or a stationary X-ray machine in a clinic, this film has what it takes to deliver a perfect film every time.

Your search is over. The 3M™ **Ultra Detail Plus** Veterinary X-ray Film is the only brand of film you need for quality, reliability, and resolution!

Double emulsion construction delivers ultimate contrast with minimum exposure



Light Crossover Problem



Depend on 3M™ Animal Care to reduce your X-ray costs

- 3M Satisfaction Guarantee
- 3M Technique guides
- 3M Technical assistance
- 3M Technician toll-free telephone support
- Cassette acquisition programs

FEATURES	ADVANTAGES	BENEFITS
<ul style="list-style-type: none"> • Double emulsion construction 	<ul style="list-style-type: none"> • Decreased exposure time 	<ul style="list-style-type: none"> • Minimizes effects of patient motion, reduces unnecessary radiation exposure, extends X-ray tube life, can be used with portable machines
<ul style="list-style-type: none"> • Zero crossover technology 	<ul style="list-style-type: none"> • Prevents light crossover, improves detail, decreased 'noise' • Provides wide range of over/under exposure tolerance 	<ul style="list-style-type: none"> • Increased diagnostic information • Fewer retakes • Reduce cost • Reduce staff & patient exposure
<ul style="list-style-type: none"> • Superior Image Detail 	<ul style="list-style-type: none"> • Reduce retakes • More diagnostic information 	<ul style="list-style-type: none"> • Saves time and money • Increased diagnostic ability
<ul style="list-style-type: none"> • Safe for all processing conditions 	<ul style="list-style-type: none"> • Convenient 	<ul style="list-style-type: none"> • No special equipment required



Animal Care Products
 3M Center, Building 270-2N-03
 St. Paul, MN 55144-1000
 1-800-848-0829

